

**GAPVD1 Rabbit Polyclonal Antibody**  
**CAB17127**



**Product Information**

<b>Product SKU:</b>	CAB17127	<b>Gene ID:</b>	26130	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse

**Additional Information**

<b>Observed MW:</b>	165kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	165kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

**Background:** Enables GTPase activating protein binding activity and guanyl-nucleotide exchange factor activity. Involved in regulation of protein transport. Located in cytosol and plasma membrane.

<b>Recommended Dilution:</b>	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
<b>Synonyms:</b>	RAP6; GAPEX5; GAPex-5; GAPVD1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GAPVD1 (NP_056450.2).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.