

# CNNM1 Rabbit Polyclonal Antibody

## CAB17131



---

### Product Information

<b>Product SKU:</b>	CAB17131	<b>Gene ID:</b>	26507	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	104kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	104kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the ancient conserved domain protein family. The encoded protein may bind copper. Alternative splicing results in multiple transcript variants.

<b>Recommended Dilution:</b>	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
<b>Synonyms:</b>	ACDP1; CLP-1; CNNM1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 692-951 of human CNNM1 (NP_065081.2).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.