

Product Information

Product SKU:	CAB17160	Gene ID:	51706	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	34kDa	Conjugate:	-
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable FAD binding activity. Predicted to be involved in bicarbonate transport. Located in extracellular exosome; membrane; and mitochondrion.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	B5R1; B5R2; B5R.1; NQO3A2; humb5R2; CYB5R1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human CYB5R1 (NP_057327.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.