

SNRK Rabbit Polyclonal Antibody



CAB17171

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

84kDa

Calculated MW:

84kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

SNRK is a member of the sucrose nonfermenting (SNF)-related kinase family of serine/threonine kinases (Kertesz et al., 2002 [PubMed 12234663]).

Immunogen information

Gene ID:

54861

Uniprot

Q9NRH2

Synonyms:

SNRK; HSNFRK


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human SNRK (NP_060189.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SNRK Polyclonal Antibody (CAB17171)