

STX17 Rabbit Polyclonal Antibody



CAB17174

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

33kDa

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55014

Uniprot

P56962

Synonyms:

STX17


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human STX17 (NP_060389.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - STX17 Rabbit pAb (CAB17174)