

FAIM Rabbit Polyclonal Antibody



CAB17176

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

20kDa/22kDa/24kDa

Applications:

IHC

Reactivity:

Human

Protein Background

The protein encoded by this gene protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene.

Immunogen information

Gene ID:

55179

Uniprot

Q9NVQ4

Synonyms:

FAIM; FAIM1

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human FAIM (NP_060617.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
