
Product Information

Product SKU:	CAB17184	Gene ID:	55966	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	-
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: Enables beta-catenin binding activity. Involved in negative regulation of cell-matrix adhesion; negative regulation of wound healing; and regulation of polarized epithelial cell differentiation. Located in several cellular components, including adherens junction; basolateral plasma membrane; and cell-cell contact zone. Is spanning component of plasma membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: MOT8; SHREW1; SHREW-1; AJAP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human AJAP1 (NP_061324.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.