

Product Information

Product SKU:	CAB17190	Gene ID:	57546	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	-
Calculated MW:	60kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene is a mitochondrial protein that functions as a phosphatase and is involved in the enzymatic resetting of the pyruvate dehydrogenase complex. Alternative splicing results in multiple transcript variants encoding the same protein.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	PPM2B; PDPC 2; PPM2C2; PDP2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 67-250 of human PDP2 (NP_065837.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.