

COL20A1 Rabbit Polyclonal Antibody



CAB17191

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

135kDa/137kDa/138kDa

Applications:

IF

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57642

Uniprot

Q9P218

Synonyms:

COL20A1; bA261N11.4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 23-290 of human COL20A1 (NP_065933.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
