

PHRF1 Rabbit Polyclonal Antibody



CAB17192

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

179kDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57661

Uniprot

Q9P1Y6

Synonyms:

PHRF1; PPP1R125; RNF221


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1350-1450 of human PHRF1 (NP_001273512.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PHRF1 Rabbit pAb (CAB17192)