

PHRF1 Rabbit Polyclonal Antibody
CAB17192



Product Information

Product SKU:	CAB17192	Gene ID:	57661	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	179kDa	Conjugate:	-
Calculated MW:	179kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable RNA polymerase binding activity. Predicted to be involved in mRNA processing and transcription by RNA polymerase II. Located in membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RNF221; PPP1R125; PHRF1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1350-1450 of human PHRF1 (NP_001273512.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.