

OSTC Rabbit Polyclonal Antibody

CAB17197



Product Information

Product SKU:	CAB17197	Gene ID:	58505	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	23kDa	Conjugate:	-
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to contribute to dolichyl-diphosphooligosaccharide-protein glycotransferase activity. Predicted to be involved in protein N-linked glycosylation via asparagine. Part of oligosaccharyltransferase complex.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	DC2; OSTC
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human OSTC (NP_067050.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.