

C12orf10 Rabbit Polyclonal Antibody



CAB17198

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

42kDa

Applications:

IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

60314

Uniprot

Q9HB07

Synonyms:

C12orf10; Gamm1; MST024; MSTP024; MYG; MYG1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 48-376 of human C12orf10 (NP_067653.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry - C12orf10 Rabbit pAb (CAB17198)