

CAB1720

Product Information

Product SKU:	CAB1720	Gene ID:	3488	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables insulin-like growth factor I binding activity. Involved in several processes, including cellular response to cAMP; regulation of smooth muscle cell migration; and regulation of smooth muscle cell proliferation. Part of insulin-like growth factor ternary complex. Biomarker of pulmonary fibrosis.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	IBP5; IGFBP5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 68-167 of human IGFBP5 (NP_000590.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.