

BRSK1 Rabbit Polyclonal Antibody



CAB17226

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

90KDa

Calculated MW:

38kDa/85kDa/86kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

84446

Uniprot

Q8TDC3

Synonyms:

BRSK1; hSAD1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 310-510 of human BRSK1 (NP_115806.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - BRSK1 Rabbit pAb (CAB17226)