

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB17247 | Gene ID: | 112464 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Rat |

Additional Information

| | | | |
|----------------|-------|------------|-----|
| Observed MW: | 35kDa | Conjugate: | - |
| Calculated MW: | 28kDa | Isotype: | IgG |

Immunogen Information

Background: The protein encoded by this gene was identified as a binding protein of the protein kinase C, delta (PRKCD). The expression of this gene in cultured cell lines is strongly induced by serum starvation. The expression of this protein was found to be down-regulated in various cancer cell lines, suggesting the possible tumor suppressor function of this protein.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | SRBC; HSRBC; PRKCDBP; cavin-3 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 75-160 of human PRKCDBP (NP_659477.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |