

C1orf131 Rabbit Polyclonal Antibody



CAB17255

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

IHC

Reactivity:

Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

128061

Uniprot

Q8NDD1

Synonyms:

C1orf131

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 13-292 of human C1orf131 (NP_001287759.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
