

PLGF Rabbit Polyclonal Antibody



CAB1727

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

25kDa/50kDa

Calculated MW:

16kDa/19kDa/24kDa/27kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a growth factor found in placenta which is homologous to vascular endothelial growth factor. Alternatively spliced transcripts encoding different isoforms have been found for this gene.

Immunogen information

Gene ID:

5228

Uniprot

P49763

Synonyms:

PGF; D12S1900; PGFL; PIGF; PLGF; PIGF-2; SHGC-10760

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 19-170 of human PLGF (NP_002623.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

