

Product Information

Product SKU:	CAB17270	Gene ID:	201626	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	67kDa	Conjugate:	-
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: Enables 3'-5'-exoribonuclease activity. Involved in several processes, including RNA metabolic process; cellular response to cytokine stimulus; and regulation of mitochondrial mRNA stability. Located in mitochondrial matrix.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	2-PDE; 2'-PDE; PDE12
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human PDE12 (NP_001309106.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.