

ATXN7L1 Rabbit Polyclonal Antibody



CAB17272

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

222255

Uniprot

Q9ULK2

Synonyms:

ATXN7L1; ATXN7L4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 450-610 of human ATXN7L1 (NP_001305158.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
