

# SLC35F1 Rabbit Polyclonal Antibody



CAB17273

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

222553

#### Uniprot

Q5T1Q4

#### Synonyms:

SLC35F1; C6orf169; dJ230I3.1

#### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC35F1 (NP\_001025029.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---