

# ZSCAN1 Rabbit Polyclonal Antibody



CAB17277

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

20kDa/45kDa

### Applications:

WB

### Reactivity:

Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

284312

#### Uniprot

Q8NBB4

#### Synonyms:

ZSCAN1; MZF-1; ZNF915

#### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ZSCAN1 (NP\_872378.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---