

ZNF530 Rabbit Polyclonal Antibody



CAB17278

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

80kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

348327

Uniprot

Q6P9A1

Synonyms:

ZNF530


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ZNF530 (NP_065931.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZNF530 Rabbit pAb (CAB17278)