

# C10orf62 Rabbit Polyclonal Antibody



CAB17284

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

IHC

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

414157

#### Uniprot

Q5T681

#### Synonyms:

C10orf62; bA548K23.1

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human C10orf62 (NP\_001009997.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry - C10orf62 Polyclonal Antibody (CAB17284)