

# RNF14 Rabbit Polyclonal Antibody



CAB17311

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9604

#### Uniprot

Q9UBS8

#### Synonyms:

RNF14; ARA54; HFB30; HRIHFB2038; TRIAD2

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 125-474 of human RNF14 (NP\_004281.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---