

MUC5A Rabbit Polyclonal Antibody



CAB17325

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

585kDa

Applications:

IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4586

Uniprot

P98088

Synonyms:

MUC5AC; MUC5; TBM; leB; mucin; mucin-5AC; MUC5A

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 5500 to the C-terminus of human MUC5A (NP_001291288.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
