

CAB17332

Product Information

Product SKU:	CAB17332	Gene ID:	4609	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human, Mouse

Additional Information

Observed MW:	57kDa	Conjugate:	Unconjugated
Calculated MW:	51kDa	Isotype:	IgG

Immunogen Information

Background:	This gene is a proto-oncogene and encodes a nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. The encoded protein forms a heterodimer with the related transcription factor MAX. This complex binds to the E box DNA consensus sequence and regulates the transcription of specific target genes. Amplification of this gene is frequently observed in numerous human cancers. Translocations involving this gene are associated with Burkitt lymphoma and multiple myeloma in human patients. There is evidence to show that translation initiates both from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site, resulting in the production of two isoforms with distinct N-termini.
Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	MRTL; MYCC; c-Myc; bHLHe39
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 347-439 of human c-Myc (NP_001341799.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.