

Product Information

Product SKU:	CAB17333	Gene ID:	9891	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	74kDa	Conjugate:	Unconjugated
Calculated MW:	74kDa	Isotype:	IgG

Immunogen Information

Background: Enables p53 binding activity and protein serine/threonine kinase activity. Involved in several processes, including protein phosphorylation; regulation of cellular senescence; and regulation of myosin-light-chain-phosphatase activity. Located in cytoplasm; microtubule cytoskeleton; and nuclear lumen.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	ARK5; NUAK1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 362-661 of human NUAK1 (NP_055655.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.