

Product Information

Product SKU:	CAB17337	Gene ID:	7827	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	-
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein that plays a role in the regulation of glomerular permeability. Mutations in this gene cause steroid-resistant nephrotic syndrome. Alternative splicing results in multiple transcript variants.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	PDCN; SRN1; NPHS2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 300-388 of human NPHS2 (NP_055440.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.