

MKI67 Rabbit Polyclonal Antibody



CAB17341

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

319kDa/358kDa

Applications:

IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

Immunogen information

Gene ID:

4288

Uniprot

P46013

Synonyms:

MKI67; KIA; MIB-; MIB-1; PPP1R105; marker of proliferation Ki-67; Ki67

Immunogen:

Recombinant protein of human MKI67

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images
