

UBA2 Rabbit Polyclonal Antibody



CAB17342

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

90kDa

Calculated MW:

60kDa/71kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10054

Uniprot

Q9UBT2

Synonyms:

UBA2; ARX; HRIHFB2115; SAE2


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 60-160 of human UBA2 (NP_005490.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images



Western blot analysis of extracts of Rat lung, using UBA2 antibody (CAB17342) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 10s.