

DHRS7 Rabbit Polyclonal Antibody



CAB17343

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

32kDa/38kDa

Applications:

WB

Reactivity:

Rat

Protein Background

This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics.

Immunogen information

Gene ID:

51635

Uniprot

Q9Y394

Synonyms:

DHRS7; CGI-86; SDR34C1; retDSR4; retSDR4

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide of human DHRS7

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images
