

# GBP5 Rabbit Polyclonal Antibody



CAB17346

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

55kDa/66kDa

### Applications:

WB

### Reactivity:

Human

## Protein Background

This gene belongs to the TRAFAC class dynamin-like GTPase superfamily. The encoded protein acts as an activator of NLRP3 inflammasome assembly and has a role in innate immunity and inflammation. Alternative splicing results in multiple transcript variants.

## Immunogen information

### Gene ID:

115362

### Uniprot

Q96PP8

### Synonyms:

GBP5; GBP-5

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant protein of human GBP5

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

## Product Images

---