

# FOXF1 Rabbit Polyclonal Antibody



CAB17356

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

40kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Protein Background

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.

## Immunogen information

### Gene ID:

2294

### Uniprot

Q12946

### Synonyms:

FOXF1; ACDMPV; FKHL5; FREAC1

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 160-379 of human FOXF1 (NP\_001442.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---