

CK1 Rabbit Polyclonal Antibody



CAB17379

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Protein Background

Immunogen information

Gene ID:

5719647

Uniprot

Applications:

Synonyms:

CK1

Reactivity:

Immunogen:

A synthetic peptide of Chlamydomonas reinhardtii CK1.

Antibody Information

Recommended dilutions:

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Product Images
