

CHCHD10 Rabbit Polyclonal Antibody



CAB17386

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

14kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a mitochondrial protein that is enriched at cristae junctions in the intermembrane space. It may play a role in cristae morphology maintenance or oxidative phosphorylation. Mutations in this gene cause frontotemporal dementia and/or amyotrophic lateral sclerosis-2. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 7 and 19.

Immunogen information

Gene ID:

400916

Uniprot

Q8WYQ3

Synonyms:

IMMD; SMAJ; MIX17A; FTDALS2; N27C7-4; C22orf16; CHCHD10

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human CHCHD10.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
