

Acetyl-FOXO1-K262/K265/K274 Rabbit Polyclonal Antibody

CAB17406



Product Information

Product SKU:	CAB17406	Gene ID:	2308	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	FKH1; FKHR; FOXO1A; Acetyl-FOXO1-K262/K265/K274
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 221-320 of human FOXO1 (NP_002006.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.