

SIAH2 Rabbit Polyclonal Antibody



CAB17423

Product Information

Size:

50uL

Observed MW:

Refer to figures

Calculated MW:

34kDa

Applications:

WB IHC IF IP

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in regulating cellular response to hypoxia.

Immunogen information

Gene ID:

6478

Uniprot

O43255

Synonyms:

SIAH2; hSiah2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50

- 1:200 IF 1:50 - 1:200 IP

1:50 - 1:200

Source:

Rabbit

Immunogen:

Recombinant Protein of human SIAH2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Isotype:

IgG

Purification:

Affinity purification

Product Images
