

CA19-9 Mouse Monoclonal Antibody



CAB17437

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

IHC IF

Reactivity:

Human

Protein Background

CA 19-9 (carbohydrate antigen 19-9, also called cancer antigen 19-9 or sialylated Lewis (a) antigen) is a tumor marker that is used primarily in the management of pancreatic cancer. CA 19-9 is an antigen defined by monoclonal antibody binding to CA 19-9, the tumor surface marker Sialyl-Lewis A.

Immunogen information

Gene ID:

Uniprot

Synonyms:

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200 IF 1:50 - 1:200

Source:

Mouse

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Human colon carcinoma cell line SW1116

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence - CA19-9 Mouse mAb (CAB17437)