

MT2A Rabbit Polyclonal Antibody



CAB17438

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

6kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4502

Uniprot

P02795

Synonyms:

MT2A; MT2


Immunogen:

Recombinant protein of human MT2A.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of HeLa cells using MT2A antibody (CAB17438) at dilution of 1:100. Blue: DAPI for nuclear staining.