

MAP1S Rabbit Polyclonal Antibody



CAB17453

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

109kDa/112kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55201

Uniprot

Q66K74

Synonyms:

MAP1S; BPY2IP1; C19orf5; MAP8; VCY2IP-1; VCY2IP1

Immunogen:

A synthetic peptide of human MAP1S.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images
