

# CDK18 Rabbit Polyclonal Antibody



CAB17515

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

54kDa

### Calculated MW:

54kDa/57kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

5129

#### Uniprot

Q07002

#### Synonyms:

CDK18; PCTAIRE; PCTAIRE3; PCTK3

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 373-474 of human CDK18 (NP\_002587.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - CDK18 Polyclonal Antibody (CAB17515)