

POU3F4 Rabbit Polyclonal Antibody



CAB17521

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

42kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the POU-III class of neural transcription factors. This family member plays a role in inner ear development. The protein is thought to be involved in the mediation of epigenetic signals which induce striatal neuron-precursor differentiation. Mutations in this gene are associated with X chromosome-linked nonsyndromic mixed deafness. [provided by RefSeq, Dec 2012]

Immunogen information

Gene ID:

5456

Uniprot

P49335

Synonyms:

BRAIN-4; BRN-4; BRN4; DFN3; DFNX2; OCT-9; OTF-9; OTF9; POU3F4


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human POU3F4 (NP_000298.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - POU3F4 Polyclonal Antibody (CAB17521)