

Product Information

Product SKU:	CAB17529	Gene ID:	6322	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable chromatin binding activity and histone binding activity. Predicted to be involved in negative regulation of transcription, DNA-templated. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SCML1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human SCML1 (NP_006737.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.