

CAB17541

Product Information

Product SKU:	CAB17541	Gene ID:	7108	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Enables delta14-sterol reductase activity. Involved in cholesterol biosynthetic process. Located in endoplasmic reticulum. Part of receptor complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ANG1; C14SR; NET47; DHCR14A; TM7SF2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 160-250 of human TM7SF2 (NP_003264.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.