

TMF1 Rabbit Polyclonal Antibody



CAB17542

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

160kDa

Calculated MW:

122kDa/123kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7110

Uniprot

P82094

Synonyms:

ARA160; TMF; TMF1


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TMF1 (NP_009045.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - TMF1 Polyclonal Antibody (CAB17542)