

## CAB17542

---

### Product Information

<b>Product SKU:</b>	CAB17542	<b>Gene ID:</b>	7110	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	160kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	123kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables androgen receptor binding activity and nuclear receptor coactivator activity. Involved in androgen receptor signaling pathway and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** TMF; ARA160; TMF1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TMF1 (NP\_009045.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.