

# ZBTB14 Rabbit Polyclonal Antibody



CAB17546

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

7541

#### Uniprot

O43829

#### Synonyms:

ZF5; ZFP-161; ZFP-5; ZFP161; ZNF478; ZBTB14

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBTB14 (NP\_003400.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---