

ZNF143 Rabbit Polyclonal Antibody



CAB17549

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

75kDa

Calculated MW:

65kDa/68kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7702

Uniprot

P52747

Synonyms:

SBF; STAF; pHZ-1; ZNF143


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 437-638 of human ZNF143 (NP_003433.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZNF143 Polyclonal Antibody (CAB17549)