

Product Information

Product SKU:	CAB17570	Gene ID:	8581	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	13kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in lymphocyte differentiation. Predicted to act upstream of or within response to stilbenoid. Predicted to be located in extracellular region and plasma membrane. Predicted to be active in cell surface.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	E48; Ly-6D; LY6D
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 25-95 of human LY6D (NP_003686.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.